sanowell

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Instructional Manual

sanowell Digi Pro 4000 Duo TENS & EMS

CE 0483

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General safety and precautionary measures. Keep the manual for subsequent use!

Do read the instruction manual, before using your sanowell Digi Pro 4000 Duo for the first time.

- Only use your TENS-EMS device according to its purpose described in the instruction manual.
- Only employ accessories recommended by hofmann gmbh.
- Persons with electronic or metallic implants (f. ex. pace maker) may not be treated by TENS or EMS without prior permission of their physician.
- In case TENS or EMS treatment is ineffective or uncomfortable stimulation has to be ceased until examination by a physician has been conducted.
- Do not apply the device when operating machinery or driving a vehicle.
- Please always turn off your TENS-EMS device before application and removal of the electrodes.
- TENS or EMS devices do not have a AP/APG-safety. They may not be operated in the environment of explosive or combustible material.
- Persons afflicted with heart disease may only apply the sanowell Digi Pro 4000 Duo under extreme caution. Please consult your physician concerning side effects before application of the device.
- TENS or EMS may not be applied during pregnancy.
- The purposeful use of TENS or EMS is nerve or muscle stimulation. If the device is applied to unsuitable body areas danger to the user may arise. Never place electrodes in the laryngeal or pharyngeal area. Muscle spasms in this area can afflict asphyxiation.
- In some rare occasions skin irritations may be caused by the electrodes. In such incidence please cease stimulation and remove the electrodes. Discontinue treatment until the cause for the skin irritation has been diagnosed.
- Do not apply electrodes to the eyes, in the mouth, or in other body holes. The electrodes are to be used only on healthy skin. Never apply electrodes to irritated or injured skin.
- Keep away from children.
- Never apply electrodes in a manner that will trans-duct currency through the brain.
- The device may not be used to treat pain of uncertain cause.
- Attention: Do not plug the electrode wires into an alternating current socket!

Additionally in Germany it is applicable as follows:

- The device is manufactured according to the medical products law (MPG).
- The device is subject to the medical products operator's regulation (MPBetreibV).
- Operator and user are obligated to strictly observe the operator's regulation in all aspects.
- For all other states of the European Community national laws and regulations for operation of medical products are valid.

Explanation of signs on the CE label

Observe the instruction manual!

Application component of type BF. Application components are not suitable to be used for direct application at the heart.

Symbol for the marking of electric appliances and electronic devices.

Manufacturer

hofmann gmbh, frankenstraße 16, 63776 mömbris
fon: +49 (0) 6029-99890, GERMANY
www.sanowell.com

Thank you very much and congratulations

We are very happy, that you chose a product from hofmann gmbh, and hope, that the sanowell Digi Pro 4000 Duo helps to relieve your disease.

Please become knowledgeable about your stimulation current unit before treatment by reading this instruction manual carefully and completely.

The sanowell Digi Pro 4000 Duo is a digital stimulation current unit with two controllable output channels for a comfortable therapy in your home. It produces electrical impulses, that may be varied through the control elements. The entire unit is built very user friendly. Nevertheless, we do recommend prior to therapy to consult your physician.

General description TENS-EMS combo unit

The sanowell Digi Pro 4000 Duo digital TENS/EMS-device is a battery-powered impulse generator that sends electrical impulses through electrodes to the body, reaching there nerves and underlying muscle groups. This device is a combination TENS unit and EMS unit, that may be used for both muscle stimulation and pain relief. The device is equipped with two independent, controllable output channels. One electrode pair can be connected to each output channel. Settings are made by means of buttons or dials.

1. General description – TENS-Function

The trans-cutic (through the skin) electrical nerve stimulation – short TENS – is a therapy method that has been used extremely successfully for decades to treat intense and chronic pain. With the sanowell Digi Pro 4000 Duo you have chosen a two-channel-TENS, that will ease or eliminate your pain through electrical impulses, that may be individually set in duration, intensity and frequency, for comfortable treatment at home.

Indication

Your physician may prescribe a TENS device for you. Please keep strictly to the doctor's orders and have him consult you in correct application of the electrodes. The correct placement is crucial to the success of your TENS therapy.

1.1. What is pain?

Under regular circumstance pain has a warning function within the human body, as they bring our attention to the fact, that something is not okay. Pain is important, as without damage to our body system would not be discovered and therefore not treated. On the other hand pain especially the chronic form puts an immense stress on the body and naturally the human soul. Nobody can and nobody should have to withstand pain for a duration of time.

Pain arises from the brain. Via sensory nerve fibres of our nervous system diversest information reach our brain to be interpreted.

1.2. What is TENS and how does TENS work?

The shortcut TENS (trans-cutic electrical nerve stimulation) depicts the stimulation of nerves by electrical impulses through the skin. In fighting against pain it offers an alternative a-thousand-fold proven in comparison with drugs, that procure a chemical intervention on the body and by that often an additional stress.

The TENS therapy bases on the natural body-own system of information transportation. Positive, as well as negative sensations (as pain) reach the brain by the same nervous system.

TENS uses the fact, that comfortable sensations travel on faster nerve strands than negative ones. That means, when the soft electrical impulses, that TENS-therapy sends out, reach the spinal cord, there is no more room left for the pain information; it can no longer be processed by the brain. Moreover TENS effects that body-own happiness hormones are spilled, the so called endorphins. These two factors together effect that the sensation of pain often is noticeably eased or completely eliminated.

A TENS-therapy may help you to get rid of your pain, so that you can face everyday business more easily. Though do consult your physician and besides pain relief by TENS occupy yourself with exploring the source.

2. General description - EMS-Function

Electrical muscle stimulation is an internationally accepted and probed treatment method for muscle injuries. It works by sending electronic impulses to the muscle to be treated, which effects passive movement of the muscles.

This is a product developed from the square pulse wave form invented by John Faraday in 1831. By use of square wave pattern it is possible to directly affect the motor sensory neurons of the muscle. EMS has low frequency, which together with the square wave pattern makes direct stimulation of the according muscle groups. It is often applied in hospitals and physical therapy clinics for treatment of muscle injuries and rebuilding of paralysed muscles in order to prevent atrophy and to improve muscle tonus and blood circulation.

2.1. How does EMS work?

EMS is applied to improve blood circulation. It loosens tight and cramped strands of muscle, stimulates muscle growth and reduces stiffness within the muscle connection. The EMS devices send comfortable impulses through the skin that stimulate the nerves within the treatment area. When the muscle receives the signal it contracts just as if it receives a signal sent from the brain. While the signal increases in strength, the muscle moves as during physical exercise. When the impulse stops the muscle relaxes and the circle starts anew (stimulation, contraction, relaxation).

3. Urinary or fecal incontinence

Incontinence is not a question of age. At least approx. 5 % of the 20- to 60-year old women and almost 15 % of the women over the age of 60 experience this bladder and bowel dysfunction. Due to the female anatomy women are considerably more often affected than men - approx. 80% of the patients suffering from incontinence are female.

Medicine differentiates the following types of incontinence:

- pressure or stress incontinence
- urge incontinence
- mixed urge- and stress incontinence
- fecal incontinence

The most common form is probably pressure or stress incontinence. This type occurs when people leak urine unintentionally even resulting from everyday pressure like coughing, sneezing, stair climbing or exercising. There are various causes for pressure or stress incontinence: nerve injuries or diseases after giving birth, infections, but also pelvic floor muscle weakness can cause stress incontinence.

The uncontrolled loss of urine, combined with partly strong to aching urinary urgency is called urge incontinence. Unlike the pressure or stress incontinence, there is mostly no sphincter dysfunction. Age- and disease-related changes in the nervous system often lead to bladder dysfunction. In many times, urge incontinence is due to inflammations or bladder outlet obstruction, such as benign prostatic hyperplasia or urethral constriction.

Fecal incontinence is a sphincter dysfunction that leads to the loss of regular control of the bowels. Involuntary excretion and leaking are common occurrences for those affected. Fecal incontinence is often caused by prostate surgery, but may also be due to longtime laxative misuse, intestinal cancer or untreated haemorrhoids.

4. What are the fundamentals of TENS-treatment?

4.1. Pulse width

Pulse width tells how long the duration of one pulse is. Through variation of intensity and pulse width it is possible to stimulate different nerve groups. So it is recommended to use different pulse width settings according to treatment mode and treated area.

You may inform yourself about pulse width settings in chapter 6.2.2.

4.2. Frequency

The frequency of the impulse in Hertz (Hz) tells in what time interval the electrical impulses are sent to the electrodes and thus to the skin. For the choice of frequency the mode of electrode placement is to be considered.

In case of stimulation in the pain area or very close to it a higher frequency (80 Hz and up), where the single impulses can not be felt anymore but treatment feels like a continuous tingling, is effective. In case of stimulation of so called trigger points (points in the muscles, that react extremely sensible) or acupuncture points a frequency of 10 Hz or less, where you can feel the single impulse, is recommended. In combination of both methods higher frequency settings are recommended. Sometimes variations of the recommended settings can be effective. You should consult your physician or try out for the most suitable setting.

You may read about setting the impulse frequency in chapter 6.2.3.

4.3. Intensity

Intensity controls the strength of the single impulse. Each patient senses the intensity individually due to tissue resistance, thickness of skin etc. Therefore intensity should be set in a range, where impulses are felt as tingling and tapping, but never as uncomfortable or even painful. Stimulation should not produce muscle contractions. These can be useful when pain is caused by cramped or spastic muscles.

Under these circumstances TENS may be used as a common muscle stimulator to

release the spasm fast. There it is recommended to set frequency high, pulse width long and intensity to visible contraction (but still tolerable for the patient). Usually 20 to 30 minutes of such muscle contraction releases the spasm. Settings for intensity can be found in chapter 6.4.1.

4.4. Treatment duration

TENS devices are commonly used for longer periods of time, 20–30 minutes. For some post-surgical treatments the duration may be prolonged to 36 hours. Pain relief partly occurs already at treatment begin. In some cases it may take up to 30 minutes until pain relieve sets in, especially when point placement of electrodes is applied or at lower frequency settings. Important is to continue pain therapy patiently and consequently.

You may read about settings for treatment duration in chapter 6.4.2.

5. What are the fundamentals of EMS-treatment?

5.1. Contraction / Relaxation

The EMS mode offers you variable muscle contraction times. At the beginning of treatment contraction times should not be set too long. A range of 1–30 sec. is available. Stimulation ensues during set time and is followed by a relaxation phase according to your individual setting.

Muscle relaxation times are also variable in this mode. A range of 1–30 sec. is available. Relaxation time should always be longer than contraction time, at least double contraction time is recommended. This cycle repeats until the timer turns off the device or you finish therapy.

You may read about setting the contraction in chapter 6.3.1. and about setting the relaxation in chapter 6.3.3.

5.2. Ascending time – Ramp

In order to prevent an uncomfortable treatment start there is the possibility to adjust the course of contractions ascending. This effects the impulse to build softly and slowly. Intensity of the electrical current reaches the highest value within the chosen time of intensity increase. In case total treatment time is shorter than the time set for intensity increase, the highest value will not be reached.



You may read about setting the ascending time in chapter 6.3.2.

5.3. Frequency

By means of the impulse frequency you may determine how many electrical impulses are sent through your skin within one second. By variation of settings the amount of running impulses per second (Hz) can be set for both channels. There are different setting options of 1 to 160 Hz.

You may read about setting the frequency in chapter 6.3.4.

6. Draft and arrangement of control elements

6.1. Control elements and their function

6.1.1. Setting into operation

For setting your sanowell Digi Pro 4000 Duo into operation turn both intensity dials clockwise until you hear a „click“ sound. By turning counter clockwise you can turn off the device.

6.2.1. Manual control setting

In case you do not want to apply one of the preset programs during TENS- or EMS-therapy, it is possible to set pulse width, frequency, and intensity manually. To do so please push first the MODE button repeatedly until the mode “N-normal treatment mode” starts blinking in the display, left. For EMS-therapy please push the MODE button repeatedly until the mode “C-continuous mode” starts blinking in the display, right. This mode allows you to set contraction time, relaxation time, ramp time and frequency manually by means of the control elements.

6.2.2. Setting the pulse width

Setting the pulse width works by pressing the WIDTH buttons. Pressing the upper button increases pulse width by 1 μ s, pressing the lower decreases it. Keeping one of the buttons pressed will increase/decrease pulse width in 5 μ s steps. You can read the pulse width setting in the centre of the display.

Recommended setting: 60 μ s

Function of pulse width is explained in chapter 4.1.

6.2.3. Setting impulse frequency

Setting of impulse frequency is achieved by pressing the RATE buttons. Pressing the upper button increases, pressing the lower decreases the frequency by steps of 1 Hz. Keeping one of the buttons pressed increases/decreases frequency in steps of 5 Hz. You can read the set frequency at the bottom of the display.

Recommended setting: 100 Hz

Function of impulse frequency can be read in chapter 4.2.

6.3. Setting in EMS-mode

For EMS-therapy please push first the MODE button repeatedly until the mode "EMS" and one of the 3 programmes "S", "C" or "A" start blinking in the display, right. Settings in the EMS-mode can be made means of the blue labelled buttons. Switching to the different settings, such as On for contraction, Ramp for ascending time, Off for relaxation and Rate for frequency, is achieved by pressing the „SET“ button. To set the several values please use the blue arrow buttons. Pressing the ENT-button activates the set values.

6.3.1. Setting contraction time - On

Setting of contraction time is achieved by pressing the „SET“ button and the blue arrow buttons. Pressing the upper button increases, pressing the lower decreases the contraction time by steps of 1 sec. You can read the set contraction time in the centre of the display.

Recommended setting: 5 seconds

Function of contraction time can be read in chapter 5.1.

6.3.2. Setting ascending time - Ramp

Setting of ascending time is achieved by pressing the „SET“ button and the blue arrow buttons. Pressing the upper button increases, pressing the lower decreases the ascending time by steps of 1 sec. You can read the set ascending time in the centre of the display.

Recommended setting: 3 seconds

Function of ascending time can be read in chapter 5.2.

6.3.3. Setting relaxation time - Off

Setting of relaxation time is achieved by pressing the „SET“ button and the blue arrow buttons. Pressing the upper button increases, pressing the lower decreases the relaxation time by steps of 1 sec. You can read the set relaxation time in the centre of the display.

Recommended setting: 12 seconds

Function of relaxation time can be read in chapter 5.1.

6.3.4. Setting frequency - Rate

Setting of frequency is achieved by pressing the „SET” button and the blue arrow buttons. Pressing the upper button increases, pressing the lower decreases the frequency by steps of 1 Hz. You can read the set frequency in the centre of the display.

Recommended setting: 100 Hz

Function of frequency can be read in chapter 4.2.

6.4.1. Setting intensity

Intensity is set by turning the intensity dials. Turning clockwise increases, turning counter clockwise decreases intensity steadily. Always start with the lowest intensity setting and work your way up slowly until you feel a tingling sensation. Treatment never ought to be uncomfortable or even painful.

The function of intensity is described in chapter 4.3.

6.4.2. Setting treatment duration

Pressing the „TIME” button increases the treatment duration by steps of 15 minutes.

15 minutes is displayed as follows: 

30 minutes is displayed as follows: 

45 minutes is displayed as follows: 

60 minutes is displayed as follows: 

You can read the set treatment duration at the bottom of the display, right.

Recommended setting: 30 minutes

Function of treatment duration can be read in chapter 4.4.

7. Preset programs

By Pressing the MODE button you can select one of seven preset TENS-programs or one of three preset EMS-programs. With each push you switch to the next lower program. The chosen TENS Program starts blinking in the display, left, the chosen EMS-Programm starts blinking in the display, right.

Setting suggestion: „M” in TENS-mode or „C” in EMS mode.

7.1. TENS-program B, or “burst treatment method”

In the first treatment method, burst (B), single bursts of 7–10 impulses individual impulses in a fixed pattern are sent. This is a combination of common TENS and low frequency TENS. In the burst mode the frequency of stimulation is set by the device and can be controlled separately through the frequency control.

7.2. TENS-program N, or “normal treatment method”

The second setting, normal treatment method (N) offers direct control of the setting elements of the device to the user: intensity, pulse width, and frequency are manually set.

How to conduct setting can be found in chapter 6.

7.3. TENS-program M, or “modulated treatment method”

The third setting is the modulated treatment mode (M). Through continuous change in intensity this method tries to prevent immunisation of the nerve. If intensity is increased in a period of low intensity within the cycle, intensity can be increased slowly by the control element, so that you will feel this as an increase

7.4. TENS-program MR, or “modulated pulse frequency”

When applying this treatment method, the pulse frequency decreases continually by 50 % and then increases until the initial value is reached.

7.5. TENS-program MRW, or “modulated pulse frequency and modulated pulse width”

Using the fifth treatment method, not only the pulse frequency, but also the pulse width are subject to increase and decrease. Oscillation occurs in the opposite manner, that means by increase of pulse frequency pulse width is decreasing, and vice versa. Intensity, pulse width, and frequency can be adjusted separately.

7.6. TENS-program SDR, or “intensity and frequency modulation”

In this sixth treatment mode intensity and frequency increase and decrease within a cycle of 6 seconds. Increase and decrease occur in an opposite manner. Intensity of stimulation decreases by 20 % while frequency increases by 45 %, before returning to the initial value. Intensity, frequency and pulse width can be adjusted separately.

7.7. TENS-program SDW, or “intensity and pulse width modulation”

The seventh treatment method contains an increase and decrease in intensity and pulse width within a cycle of 6 seconds. Increase and decrease occur in opposite oscillation. Stimulation intensity decreases by 20 % while pulse width increases by 54 % until they return to the initial value. Intensity, pulse width, and frequency are separately adjustable.

7.8. EMS-program S, or “Dual-Mode”

Stimulation on both channels is synchronal. “Contraction time” includes „ramp up” and „ramp down” time. When setting the “contraction time”, please therefore observe that “contraction time” should be at least double „ramp”-time in this mode.

7.9. EMS-program C, or “Continuous Mode”

The second setting, normal treatment method (C) offers direct control of the setting elements of the device to the user: contraction time, relaxation time, ramp time and frequency are manually set.

7.10. EMS-program A, or „Alternating Mode”

Stimulation on channel 2 starts after the first contraction on channel 1 has finished. When setting the “contraction time”, please observe that “contraction time” should be at least double „ramp”-time in this mode. “Relaxation-time” should be the same or more than “contraction-time”.

8. To connect a vaginal- or rectal probe to the sanowell Digi Pro 4000 Duo

- Please read the instruction manual of the used vaginal- or rectal probe as well as the instruction manual of the hand control carefully and completely before treatment.
- Connect the electrode cable to vaginal- or rectal probe.
- Push the plug of the Digi Pro 4000 Duo hand control into the socket marked „Probe” on the side of the device and activate the hand control mode at the device. You may read about how to make settings in chapter 8.2.
- Spread a light coating of water based lubricating gel over the main body of the probe and insert it into the corresponding body hole (water based lubricating gel can be obtained from suppliers of medical products or pharmacies).

8.1. To use the hand control

The hand control mode offers you additional ease of use for vaginal- or rectal probes. After you have made all the settings at the device, you can switch on or switch off the current flow through the probe easily and conveniently by pushing the button at the hand control. By pressing the red button at the hand control current flows through the probe or electrode, releasing the red button stops the current flow. This quick and direct possibility for comfortable operation provides more security and additional comfort to the patient.

8.2. To set the hand control mode at the device

To access the hand control mode please press the „Time“-button and keep it pressed while you set the device into operation by turning one intensity dial clockwise.

8.3. To activate the hand control mode at the device

To activate the hand control mode press the „Mode“ button first, and then the „Rate“ (arrow up) button. For checking purposes, an „M“ appears in the top left-hand corner of the display.

8.4. To deactivate the hand control mode at the device

Press the „Mode“ button first, and then the „Rate“ (arrow up) button. The „M“ in the top left-hand corner of the display disappears.

8.5. To restart the device

Turn the device off for a short time by means of the intensity dial and then turn it on again. Now, the device is ready to operate again.

8.6. Recommended settings to be used with a vaginal- or rectal probe:

Device mode: EMS

Contraction: approx. 3 – 5 sec.

Ascending time: at the beginning of treatment 3 sec.

Relaxation: approx. 10 – 12 sec.

Frequency: 100 Hz

9. What to observe about the electrodes

Your physician ought to decide, which electrodes to best use for your kind of pain. Do observe the instructions on electrode placement in the instruction manual in order to keep stimulation effective during treatment time and to prevent skin irritations. You may only use electrodes with the CE-certification sign.

The wires delivered with the TENS device are plugged into the sockets at the device. To do so place the device in your palm and push the plugs at the end of the wire into one of the sockets (see drawing, chapter 6); you may use one or two connection wires. After connecting the wires to the stimulation unit please connect each electrode to one wire end (positive and negative terminal). Please be careful when attaching or disconnecting the wires. Do not pull at the wires, as this may lead to breakage.

9.1. Electrode placement

Best possible electrode placement is the key to successful treatment. Besides the most usual placement strategies described as follows, it may be effective to use a deflective electrode application. Thus it is important that you take your time in experimenting for the best electrode placement for you, as the positioning is individual. Please have a TENS firm physician consult you.

9.1.1. The bordering placement

This is the most common placement strategy. Here the painful area is bordered by four electrodes. You can decide whether the current is to flow directly through the pain area or surrounds it. In order to obtain a direct flow through the spot of pain you have to place the positive terminals of the electrode pairs diagonally to the negative ones; if trying to achieve a surrounding flow, they have to be placed parallel. This is very effective if the pain is in an extremity or deep within the tissue.

9.1.2. Placement in the area of a spinal nerve

The spinal nerves exit the spinal cord in pairs between the vertebrae. In this technique of electrode placement stimulation happens above the innervated area by placing one electrode exactly on the spot of pain and the other at the nerve's root at the spinal cord.

9.1.3. Placement at motive, trigger, and acupuncture points

Motive, trigger, and acupuncture points are places of high tissue conductivity. In order to destine the exact position of such a point please consult your physician.

9.1.4. Multiple placement strategies

As your TENS device possesses two independent channels, you can apply different strategies in electrode placement at the same time. One channel can be used for direct stimulation of the pain area, whereas the second channel is applied for point therapy by bordering placement.

9.2. Use of self-adhesive electrodes

9.2.1. Application

- Clean the skin thoroughly with soap and water before application of the electrodes
- Excessive hair should be cut, but not shaved
- Connect the wires to the electrodes
- Pull the safety cover off the electrodes and apply firmly to the skin area by pressing.

9.2.2. Removal

- Lift the electrodes at the corner and gently pull off. Do not pull at the electrode wire, as this might damage the electrodes
- Re-glue the electrodes to the safety cover and then remove the wires by carefully turning and pulling.

9.2.3. Care and storage

Your sanowell Digi Pro 4000 Duo will only achieve best results when electrodes and wires are cared for accordingly.

- Clean the wires with a damp cloth; slight powdering with talcum (f. ex. baby powder) will prevent tangling of the wires and improves their durability
- After each application please store the electrodes in their re-closeable bags in a cool dry place (f. ex. the refrigerator)

- Life time of the electrodes is improved by applying a drop of cool tap water (please do not use distilled water) to the adhesive area and air drying the surface turned up. Moistening the adhesive area too much may spoil the adhesive capability.

10. Helpful hints for skin care

In order to prevent skin irritations we suggest as follows:

- Wash the areas where the electrodes are to be placed prior to application and after removal with a mild soap and water. Rinse the soap thoroughly and dry the skin.
- Excessive hair growth may be cut away with the scissors; do not shave the area that is to be stimulated.
- Many skin problems arise due to the electrodes being pulled too tight over the skin at application. To prevent thus apply the electrodes from their centre to the rim and do not pull them over the skin.
- To prevent irritation by stretching tape the excess wires in a circle to the skin so it will not pull on the electrodes.
- When removing the electrodes always do so in the direction of hair growth.
- It may be helpful to apply skin lotion on the areas that electrodes were attached to between applications.
- Never adhere electrodes to skin already irritated or injured.

Attention

- Do not apply to injured skin.
 - Dispose of electrodes that do not adhere well.
 - Self-adhesive electrodes are to be used by the same person only.
 - In case of skin irritation cease usage and consult your physician.
 - Prior to application read the instruction manual for self-adhesive electrodes.
- If stimulation is uncomfortable either reduce intensity to a tolerable amount or discontinue stimulation and ask your physician.

11. Handling the Digi Pro 4000 Duo

Your sanowell Digi Pro 4000 Duo is already equipped with a battery, but may also be operated by storage batteries.

11.1. Usage of re-chargeable batteries

Prior to usage of a new device the storage battery has to be charged according to the manufacturer's prescription. Please read all instructions and safety precautions printed on the accumulator and in this leaflet before usage. Stored for more than 60 days the re-chargeable battery may discharge. Therefore prior to usage they need to be charged again, if they were stored for a duration of time.

11.1.1. Charging the accumulator

- Plug the accumulator into a 220-230V socket. Using an accumulator not delivered as accessory may cause fire, electric shock or body injury.
- Regarding charging time observe the instructions of the manufacturer.
- After elapse of charging time unplug the accumulator and remove the charged battery.
- Re-chargeable batteries should always be stored loaded.

To guarantee best battery performance the following rules should be observed:

- Overloading of up to 24 hours does not damage the re-chargeable battery. Though frequent overloading may reduce their life time.
- Do not short plug the accumulator terminals. This causes the accumulator to run hot and may cause permanent damage. Do not keep batteries in clothing or purses as terminals may contact coins, keys, or other metal objects.

Safety precautions:

- Do not charge any other batteries with the accumulator than re-chargeable batteries. Other tape batteries may leak or explode.
- Do not burn re-chargeable batteries, as they may explode!

11.2. Testing and exchanging batteries

Delivery already contains a 9V battery. If the „low battery“  sign shows in the LCD display, the battery is too weak to operate the device and needs to be exchanged. The device automatically turns off and can not be turned on again until a new battery has been inserted. In order to guarantee save device function from time to time the battery needs to be replaced.

When replacing the battery proceed as follows:

- Pull off the battery compartment's lid downwards.
- Remove the consumed battery from the battery compartment.
- Insert a new 9V battery according to the drawing in the bottom of the battery compartment. The destined polarity has to match.
- Push the lid back upwards.

11.3. Care, transportation, and storage

- The device may be cleaned with alcohol. Attention: When handling combustible fluids do not smoke or operate with open flame (e.g. candles, etc.)
- Stains and massive soiling may be removed with a mild detergent.
- Do not immerse the device into liquids or subject it to large amounts of water.
- For safe transportation always insert device into the carrying case.
- If the device is not used for a prolonged period of time, remove the batteries from the battery compartment (acid could leak and damage the device). Store the device and all accessories cool and dry in the carrying case.

11.4. Safety-technical control measures, STC

According to the medical products operator's regulation (MPBetreibV), § 6 (safety-technical control measures), the user is obliged to have regular safety-technical controls conducted. hofmann gmbh prescribes these controls in accordance to MPBetreibV § 6 in a 24-months turn. The safety-technical control may only be delegated to persons, who are able to conduct controls in an orderly fashion based on their education, knowledge, and experience achieved by practical work, and who also have suitable measuring and proofing tools available.

Standards for safety-technical control are to be requested from the quality management agent of hofmann gmbh.

Prior to each application of the sanowell Digi Pro 4000 Duo device the following controls are to be conducted:

1. Inspect device for surface damage:
 - Deformation of casing?
 - Damaged or defect electrode sockets?
2. Inspect device for damaged control elements:
 - Imprints still legible?
3. Inspect display elements:
 - The LCD-display has to show device settings when turned on.
4. Inspect accessories for usability:
 - Electrode wires non-damaged?
 - Electrodes non-damaged (no fissures, no worn material)?

11.5. Malfunctions

Should malfunctions occur, when using the Digi Pro 4000 Duo, please check as follows:

- Is the treatment mode set correctly? If not, please correct setting.
- Is the wire connected correctly to the device? Wires are to be completely inserted into the sockets.
- Does LCD-display show the device settings? If necessary, exchange the battery.
- Is the wire damaged? In any case of damage, exchange the wire.
- Are 2 Electrodes connected to the paired wire? In case of any other problems, return the device to the dealer. Do not try to repair a defect device yourself! Manipulation voids the warranty!

For problems with the device or accessories please contact your dealer.

11.6. Disposal

Device packaging should be disposed of in the recycling bin, metal parts of the device in the metal recycling disposal. Plastic parts, electrical parts, and electronic chips are to be disposed of as electronic garbage. Batteries are problem garbage and are to be disposed of accordingly. Return consumed batteries to your local recycling plant. Do not dispose of in regular garbage. You may obtain further information on disposal of problem material at your local waste management centre.

12. Contents of delivery and accessories

12.1. Contents of delivery

- 1 x sanowell Digi Pro 4000 Duo TENS & EMS unit
- 2 x lead wires
- 4 x self-adhesive electrodes, size 40 x 40 mm
- 1 x battery, 9 V
- 1 x instruction manual
- 1 x carrying case

Caution:

The above mentioned contents are standard equipment, that can vary within distributors.

12.2. Accessories for incontinence therapy

- hand control for the sanowell Digi Pro 4000 Duo, article No. 89243
- sanowell vaginal probe, for the treatment of urinary incontinence in women, article No. 89241
- sanowell rectal probe, for the treatment of urinary incontinence in men or fecal incontinence in women and men, article No. 89242

12.3. Accessories

- self-adhesive electrodes, 2mm plug, 40 x 40 mm, article No. 89230
- self-adhesive electrodes, 2 mm plug, 40 x 80 mm, article No. 89231
- self-adhesive electrodes, 2 mm plug, size 32 mm round, article No. 89232
- permanent electrodes for use with electrode gel, article No. 89233
- electrode gel for TENS and EMS devices, article No. 89240
- cables for permanent electrodes, article No. 89245
- cables for self-adhesive electrodes, article No. 89244

If you would like to make your device compatible with the permanent electrodes, you need to supply yourself with the fitting cables and the electrode gel.

13. Technical specifications

Channels:	dual, isolated between channels
Pulse Amplitude:	adjustable to maximal 80 mA for positive impulse or to about 10 mA peak for negative impulse at 500 Ohm load per channel.
Pulse Rate:	adjustable from 1 to 160 Hz in steps of 1 Hz
Pulse Width:	adjustable from 50 to 260 μ s in steps of 1 μ s.
Contraction time:	1 sec. - 30 sec. (adjustable)
Relaxation time:	1 sec. - 30. sec (adjustable)
Ramp time:	1 sec. - 8 sec. (adjustable)
Timer:	15, 30, 45 and 60 minutes
Dual flank control by software:	increase of amplitude linear at change of modus (anti shock)
Treatment duration counter:	shows time of device usage by patient
Keylock:	prevents manipulation of parameters set by physician through user
Wave Form:	asymmetrical bi-phasic square pulse
LCD-display:	zur Anzeige der eingestellten Parameter
Power Source:	9V block battery
Operating Conditions:	approximately 60 hours at regular settings
Betriebsbedingungen:	0 °C to +50 °C at 20 to 75% relative humidity, however with no condensation, pressure 700 hPa to 1060 hPa.
Storage Conditions:	0 °C to +50 °C at 20 to 75% relative humidity, however with no condensation, pressure 500 hPa to 1060 hPa.
Dimensions:	128 mm (height) x 70 mm (width) x 25 mm (length)
Weight:	170 g (inclusive of battery)

Safety for electrical shock, type BF.

All electrical specifications valid at ± 20 % 500 Ohm load.

14. Warranty

We give a 2-year warranty for our products, starting with the day of purchase.

Within this warranty period, we will eliminate defects caused by material or production faults free of charge. We will eliminate defects either by repair or by replacing individual parts or the entire device.

The warranty does not cover defects that are caused by improper use of the device (e.g. connection to the wrong type of current or voltage, connection to unsuitable power sources, breakage due to the device falling down etc.). Normal wear and defects that represent a negligible impairment of the device's value or serviceability are also not covered by the warranty.

The warranty expires in case of unauthorized manipulation or use of other than original spare parts.

In case of a defect that is covered by the warranty, please send the device and the original receipt to one of our service partners or directly to hofmann gmbh.

Any other claims based on this warranty are excluded unless we are legally obliged to accept liability. Warranty claims against the seller resulting from the sales contract are unaffected by this warranty.

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